

Available Communication Devices



UBIDUO

Purpose: Text with someone in real-time

Quick Use Steps:

- Turn on both units
- Type and communicate in real time



POCKETALKER

Purpose: Use to boost important sounds and reduce background noise

Quick Use:

- Plug in the headset
- Turn on the device and adjust volume




PORTABLE HEARING LOOP

Purpose: Use to send sounds directly to your hearing aids and cochlear implant

Quick Use:

- Connect device to hearing aid or cochlear implant
- Ensure hearing aid/cochlear is set to T-Coil mode



GOTALK

Purpose: Express yourself using recorded messages for symbols and icons

Quick Use:

- Hold down record button to record your message
- Select the desired word and/or picture icon(s)
- Press to play recorded message



VIDEO PHONE (VP) FOR DEAF USERS

Purpose: Use to make a video call to hearing individuals

Quick Use:

- The video call connects to a qualified ASL interpreter
- Interpretation is provided between sign (for the deaf person) and spoken language (for the hearing caller)

*No cost to the caller